


ARTICULATED ARM

Articulated arm for arthroscopic and open surgery
shoulder, knee, elbow, ankle and small joints



DESCRIPTION

The removable articulated arm can be fixed to the operating table whatever the patient's position: lateral decubitus, dorsal decubitus, beach chair position, etc. It provides effective and stable support in all positions to ensure surgical comfort.

Stable fixation to the operating table

Rotation **360°**

Quick and **easy connection** of supports

Control handle

Simple handling without electrical or pneumatic energy



ACCESSORIES

A wide variety of highly versatile accessories are available for use with the removable articulated arm , to adapt quickly and easily for many type of surgical procedures.



Hand support



Connection Part



Traction support



Foot support




Leg support



Straight arm support



Vertical arm support



Hip support



Traction strap



Storage box

DESCRIPTION	REFERENCE
Articulated arm	VSPDR-0001
Connection Part	VSPDR-0001-01
Straight arm support	VSPDR-0001-02
Foot support	VSPDR-0001-03
Vertical arm support	VSPDR-0001-04
Hip support	VSPDR-0001-05
Leg support	VSPDR-0001-06
Hand support	VSPDR-0001-07
Traction support	VSPDR-0001-08
Traction strap	VSPDK-0001-08
Storage box	VSPDR-0000



