




# Zimmer Biomet Equipment Cart

Brochure





\*varies depending on cart configuration



	<b>Cart 27-inch PRO</b>	<b>Cart 32-inch PRO</b>	<b>CART 2x27-inch PRO</b>
<b>Item Code</b>	11010036	11010037	11010038
<b>Size and Weight</b>	Height: 1466 mm Width: 618 mm Depth: 624 mm Total Weight: 116 kg	Height: 1466 mm Width: 618 mm Depth: 624 mm Total Weight: 128 kg	Height: 1466 mm Width: 618 mm Depth: 624 mm Total Weight: 128 kg
<b>Maximum Load</b>	Total: 180 kg, Shelves: 50 kg, Drawer: 3kg		
<b>Shelves &amp; Drawers</b>	4 x Std. Shelf 1 x Pull-Out Shelf with Handle 1 x Drawer 1 x Keyboard Drawer	4 x Std. Shelf 1 x Pull-Out Shelf with Handle 1 x Drawer 1 x Keyboard Drawer	4 x Std. Shelf 1 x Pull-Out Shelf with Handle 1 x Drawer 1 x Keyboard Drawer
<b>Dimensions Shelves &amp; Drawers</b>	Shelf: 450 x 491 mm, Total Load: 50 kg Drawers: 450 x 491 x 96 mm, Total Load: 3 kg		
<b>Monitor Holder</b>	VESA 75/100 Tilt of +/- 15° Swivel of +/- 35° Max. Load: 0 – 14 kg	not available	VESA 75/100 Tilt of +/- 15° Swivel of +/- 80° Max. Load: 0 – 14 kg
<b>Variable Height Arm</b>	not available	VESA 75/100 Tilt of +/- 20° Rotation of +/- 5° Swivel of +/- 80° Max. Load: 11–18 kg Length: 300 mm	VESA 75/100 Tilt of +/- 20° Rotation of +/- 5° Swivel of +/- 80° Max. Load: 8–14 kg Length: 300 mm
<b>IV Pole</b>	4-arm safety bottle holder, Load capacity: 4 x 2 kg Height-adjustable range: 870 mm		
<b>Camera Holder</b>	Yes	Yes	Yes
<b>Cable Winder</b>	Yes	Yes	Yes
<b>Equipment Fixing</b>	Yes	Yes	Yes
<b>Isolating Transformer</b>	Transformer power: 2000 VA, Line voltage: 230 V, Frequency: 50/60 Hz, without insulation monitor*		
<b>Equipotential Bonding</b>	Panel with 8 equipotential plugs		
<b>Castors</b>	4 twin swivel castors, lockable, Ø 125 mm, conductive		
<b>Cable Configuration</b>	Equipped with 1 x power cable EU (Schuko), 5m and 5 x IEC-to-IEC Euro cables		Equipped with 1 x power cable EU (Schuko), 5m and 6 x IEC-to-IEC Euro cables

\* Insulation monitor is available as a separate item code

## Other Accessories

Item Code	Description
11010083	Insulation Monitor Upgrade Kit
11010084	Gas Bottle Holder, Dia. 145mm, max. 20 Kg
11010085	Bottle Holder, Dia. 120mm, max. 5 Kg
11010086	Bottle Holder, Dia. 145mm, max. 5 Kg
11010087	Carrier Arm VESA 75/100 for Tablet Holder
11010089	Mounting Kit for Carrier Arm (REF 11010087)

## Power Cables and Equipotential Cables

Item Code	Description
11010039	Power Cable UK with equipotential bonding cord, 5 m
11010090	Power Cable UK, 5 m
11010040	Power Cable EU with equipotential bonding cord, 5 m
11010041	Power Cable CH with equipotential bonding cord, 5 m
11010058	Equipotential Cable, 0.3 m
11010059	Equipotential Cable, 0.5 m
11010060	Equipotential Cable, 1.0 m
11010061	Equipotential Cable, 1.5 m
11010094	Equipotential Cable, 3.0 m
11010053	IEC-to-IEC Euro Cable, 0.6 m
11010054	IEC-to-IEC Euro Cable, 1.0 m
11010055	IEC-to-IEC Euro Cable, 1.5 m
11010056	IEC-to-IEC Euro Cable, 2.0 m
11010057	IEC-to-IEC Euro Cable, 2.5 m

**For more information please contact your Zimmer Biomet representative.**

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet. This material is intended for health care professionals. Distribution to any other recipient is prohibited. Not for distribution in France. For product information, including indications, contraindications, warnings, precautions, potential adverse effects and patient counselling information, see the instructions for use and [www.zimmerbiomet.com](http://www.zimmerbiomet.com).

©2018 Zimmer Biomet



**ZIMMER BIOMET**  
Your progress. Our promise.®

1085.2-EMEA-en-REV0118

Distributed by  
Zimmer GmbH  
Sulzerallee 8  
8404 Winterthur  
Switzerland

[www.zimmerbiomet.com](http://www.zimmerbiomet.com)



Legal Manufacturer  
ITD GmbH  
Sportplatzstraße 3  
84381 Johanniskirchen  
Germany

[www.itd-cart.com](http://www.itd-cart.com)

