


# EZPass Suture Passers

Braided Nitinol has been optimized and combined with an oxide coating to minimize the friction transferred to soft tissue.

Seven different curvatures to provide multiple options from any portal position.

Sharp, low profile tips to enhance the ease of tissue penetration while reducing tissue trauma.



Thumb wheel allowing for quick and easy advancement of the Nitinol passing kite.

Disposable instrument packaged sterile and ready for surgical use.

An ergonomic handle design allows the use of a natural, relaxed, grip.

Unique colors allowing for clear identification of desired style.

# BIOMET

## Ordering Information

Part Number	Description
110008593	90° Up 10pk
110008590	45° Up 10pk
110008589	45° Dog Leg Up 10pk
110008594	70° Right 10pk
110008588	70° Left 10pk
110008591	30° Right 10pk
110008592	30° Left 10pk

Part Number	Description
110007380	Braided Nitinol Kite
904050	90° Up
904051	45° Up
904052	45° Dog Leg Up
904053	70° Right
904054	70° Left
904055	30° Right
904056	30° Left



This material is intended for health care professionals and the Biomet sales force only. Distribution to any other recipient is prohibited. All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to Biomet Inc. or its affiliates unless otherwise indicated. This material must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Biomet.

Check for country product clearances and reference product specific instructions for use. For complete product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see the package insert, [www.Biomet.com](http://www.Biomet.com), or contact your local Biomet representative.

Not for distribution in France.

**BIOMET**  
SPORTS MEDICINE  
One Surgeon. One Patient.

**Legal Manufacturer**  
Biomet Sports Medicine  
56 East Bell Drive  
P.O. Box 587  
Warsaw, Indiana 46581  
USA