

A **stronger**¹ anchor with a *softer*² touch

Introducing the

JuggerStitch[™] All-suture, Knotless Meniscal Repair Device

- Ergonomic delivery mechanism that allows for **quick and simple** implant deployment
- The low profile **taper tip needle** intended to minimize cutting the meniscal tissue during insertion
- The **all suture knotless implant** allows for controlled tensioning while eliminating rigid plastic anchors and hard meniscal surface knots

For more information,
ask your Zimmer Biomet representative
or visit us at zimmerbiomet.com



1. *Static Loading of JuggerStitch Meniscal Repair Device*. BSM26.VR2. 24 AUG 2016.
2. Compared to current hard PEEK implants.

