


JuggerKnot™ SOFT ANCHOR

The **JuggerKnot™ Soft Anchor—1.0 mm Mini** represents the next generation of suture anchor technology. The 1.0 mm deployable anchor is completely suture-based and the first of its kind. With over 75,000 anchors sold, the award-winning JuggerKnot™ Soft Anchor is now designed specifically for hand and wrist procedures. The new configuration includes two T-28 tapered needles to aid in soft tissue reattachment.

The **JuggerKnot™ Soft Anchor—1.0 mm Mini** is indicated for:



- Repair/Reconstruction of collateral ligaments
- Scapholunate ligament reconstruction
- Ulnar or lateral collateral ligament reconstruction
- Flexor and extensor tendon at the PIP, DIP, and MCP joints for all digits.



100% suture anchor

Available with either
#2-0 or #3-0 suture

#1 MaxBraid™ Suture
Polyester Sleeve



26 lbs of pull-out strength for #2-0
suture and 20 lbs of pull-out strength
for #3-0 suture¹

Sleeve protector
eliminates need for
drill guide

Two T-28 tapered needles

Drill hole is
only 1 mm

1. Data on file at Biomet, Inc. Bench testing results are not indicative of clinical performance.