

**CLASSIC TOWER** 



THE MOST ADVANCED VISUALIZATION SYSTEM









# OUR GOLD STANDARD SYSTEM, COMPACT & DESIGN






OpenVcam









## **HIGH IMAGE QUALITY 4K**

- 4K UHD Resolution 2160p
- 32" IPS screen

## **ALL-IN-ONE PROCESSOR**

- Mono-led video light source (60 000h)
- Recording picture & video on a Flashdrive 3.0
- Multi-display: until 3 displays at the same time
- Up to 3 videos outputs (Display Port)
- Automatic settings & configuration of white-balance

## **BROADCAST & VISIOBROADCAST**

- Recording of surgeries on a dedicated server
- Retransmission of surgeries on a webpage for the surgeon
- Broadcast of live video & audio (bilateral sound)

#### **RELATED CAMERAS**

- 1. Monobloc cameras Associated with a sterile single-use cover.
- Directional system with a rotation angle of 170° to -170° for 30° cameras
- Ergonomic handle with Photos/videos & light modes buttons
- A single cable for video & light sources
- ArthroVims 4mm 30°: camera for arthroscopy
- OpenVcam on articulated arm : environment camera to film open surgeries

#### 2. Proximal camera

- EndoVcam
- Compatible with all rigid scopes under Ø5mm
- 4K UHD image quality.

# **ENDOVCOVER SINGLE-USE STERILE COVER**

- Very resistant, hermetic
- Ultra-fast fitting (50s)
- No condensation & contamination risks
- Maintain high quality images
- Sterilization free for light cable and scope







REFERENCE	DESCRIPTION
	CLASSIC CONFIGURATION
VLS750-4H	Video light processor
PVM594	32" IPS screen
CTH0,	Medical Vcart
	CAMERAS
VSX4030J-4F	ArthroVims Ø4 mm - 30°
VSX8000PR-4I	EndoVcam - Proximal camera
VSX8030S-4I	OpenVcam with articulated arm
	STERILE SINGLE-USE COVER
LW4030J-UHI	EndoVcover for VSX4030J-4K
LW10030-S-UHI	EndoVcover for VSX8030S-4K
	OPTIONS
OPN0	Vcart option with articulated arm
DH3:	Vcart option for dual monitor 32"
BR0	Broadcast with Touchless interface
BR0: CTH05-3:	Broadcast Veart for additional monitor
UR-4MI	External recorder with Dicom
On-4ivit	External recorder with bloom
CLX9000-:	ACCESSORIES
CLX9000-:	Proximal camera straight light cable Proximal camera angled light cable
10001409-0	Trocar sheath Shoulder for ArthroVims
10001409-0.	Trocar sheath Knee for ArthroVims
1000140	Trocar obturator for ArthroVims
04-20-02	USB Keyboard



Articulated Arm



Dual Monitor



Additional Screen

